

Further information at
www.gpescientific.co.uk















BVC-accessories - at a glance

Accessories for fluid aspiration systems
BVC basic, BVC control and BVC professional

Accessories for fluid aspiration systems

Biochem-VacuuCenter

The compact BioChem-VacuuCenters offers perfect solutions for safe, sensitive and efficient aspiration of supernatants in all application areas. All systems are available either with vacuum-resistant, autoclavable 4L PP or a 2L borosilicate glass bottles, with glass preferable for working with aggressive disinfectants such as chlorine bleach. All versions are ready for connection of a second hand control. A range of accessories are available:

| | |
|--|--|
| <p>Handset VHC^{pro} NEW</p> <p>Ergonomic light-weight aspiration handset that comes with a 2.5 m flow through, a very flexible silicone tube. Fully autoclavable.</p>  <p>688061</p> | <p>BVC shuttle</p> <p>Mobile stand for BVC-models (Control and Professional). Designed for easier cleaning and changing.</p>  <p>696880</p> |
| <p>Table stand VHC^{pro} NEW</p> <p>Stable stand for the VHC^{pro}</p>  <p>635680</p> | <p>4L bottle (PP)</p> <p>Chemically resistant 4L collection bottle. Designed with a special bottle cap that reduces aerosol and foam formation.</p>  <p>635810</p> |
| <p>Tip ejector VHC^{pro} NEW</p> <p>Adapter with ejector for pipette tips 2 – 200 micron.</p>  <p>635638</p> | <p>2L bottle glass</p> <p>2L borosilicate glass bottle. Designed with a plastic coating in order to prevent biohazardous situations.</p>  <p>635809</p> |
| <p>Extension kit second VHC/VHC^{pro} connection</p> <p>Hose nozzle and inlet tube for minimised formation of aerosol and foam (VHC/VHC^{pro} not included).</p>  <p>699943</p> | <p>Adapter for 2L glass bottle</p> <p>Adapter for 2L-glass bottle (BVC Control, BVC Professional).</p>  <p>635839</p> |
| <p>Quick couplings VHC/VHC^{pro} - bottle</p> <p>PVDF quick couplings with adapter for VHC/VHC^{pro} connection to bottle. Also usable as an extension kit (instead of using 699943).</p>  <p>635807</p> | <p>Spare tubing for VHC/VHC^{pro}</p> <p>Silicone tubing (per metre).</p>  <p>636156</p> |
| <p>Quick coupling set connection bottle-pump</p> <p>PVDF quick couplings for connection of bottle to pump BVC or BVC basic.</p>  <p>635808</p> | <p>Protection filter 0.2 micron</p> <p>Hydrophobic protection filter is designed to protect the pump, electronics and the environment. 20-times autoclavable.</p>  <p>638266</p> |

GPE SCIENTIFIC LTD

Greaves Way Industrial Estate, Stanbridge Road, Leighton Buzzard, Bedfordshire, LU7 4UB, United Kingdom
 T +44 (0)1525 382277 · F +44 (0)1525 382263
 info@gpescientific.co.uk · www.gpescientific.co.uk