


# SNOWBIRD -90°C FREEZER

## SPECIFICATIONS

Supplied by



Perfect for:

- Small volume samples
- Isolation of sensitive or infectious samples
- Irreplaceable or one-off materials
- Samples that require high-security storage

For more information please visit:

[www.gpescientific.co.uk](http://www.gpescientific.co.uk)

SnowBird	Unit	Part no: 8005121010
Temperature range	°C	-65 to -90
Temperature stability	°C	+/- 2,5
Noise level	dB(A)	< 52
Volume internal	L	87
Weight	Kg	106
External dimensions (DxWxH) - With castors - Adjustable feet mounted	mm	716x650x863
External dimensions (DxWxH) - Without castors - Adjustable feet mounted	mm	716x650x850,5
Internal dimensions (DxWxH)	mm	510x440x390
Shelf dimensions (DxW)	mm	475x440
Cooling system		Cascade
Refrigerant 1. stage		R507
Refrigerant 2. stage		R170
Insulation	mm	100
Power consumption average	W	400
Power consumption maximum	W	980
Voltage	V/Hz	230/50

GPE Scientific Limited, Greaves Way Industrial Estate, Stanbridge Road,  
Leighton Buzzard, Bedfordshire, LU7 4UB, United Kingdom

T: +44 (0) 1525 382277  
F: +44 (0) 1525 382263

sales@gpescientific.co.uk  
www.gpescientific.co.uk



Manufactured by

