

# SCHLENK LINE COMPLETE SET UP

We are offering a package deal on GPE Schlenk Lines purchased with a Vacuubrand VacuuView extended Vacuum Gauge (683210) and a Vacuubrand RZ 2.5 Rotary Vane Pump (698122).

We can customise the Schlenk line to suit your requirements with either J. Young Taps or High Vacuum Glass Taps, Hose Barbs or Ground Glass Joint Connections.

## EXCLUSIVE DEAL - SAVE OVER £300

on a Vacuubrand Pump and Gauge when purchasing any GPE Schlenk Line

QUOTE - SCHLENKVACDEAL when placing an order



698122 - Vacuubrand RZ 2.5 Rotary Vane Pump  
Two stage, small flange with  
hose nuzzle DN 10 mm for inlet  
230V /50-60 Hz  
UK mains Cable




683210 - Vacuubrand  
Vacuuview extended Gauge  
110-0.001 mbar  
100-230v/50-60Hz  
UK Plug

Pricing on Vacuubrand items listed above is only valid when purchasing a GPE Scientific Schlenk line.  
Discount is not available on items other than those listed above.

Contact Us

## SCHLENK LINE COMPLETE SET UP



### SCHLENK LINES

*Our range of Schlenk Lines is exclusive to GPE Scientific.*

*We manufacture our Vacuum Manifolds / Schlenk Lines using the official J. Young High Vacuum PTFE valves alternatively these can be manufactured with glass stopcocks.*

*We produce our standard range of Schlenk Lines which includes General Purpose, Space Saving and Dual Bank Manifolds.*

*We can also custom make Schlenk Lines to suit your needs.*



### RZ 2.5 ROTARY VANE PUMP

*High water vapour tolerance due to efficient gas ballast; very good ultimate vacuum even with gas ballast.*

*Vacuum-tight at switch-off; external anti-suckback valve not needed.*

*Large oil volume: Long intervals between oil changes and ease of maintenance due to telescopic design.*



### VACUUM GAUGE

*VACUU-VIEW Extended Gauge*

*Precision and chemical resistance in an exceptionally wide range from atmosphere down to 10<sup>-3</sup> mbar, one gauge only for both rough and fine vacuum.*

*The heavy duty combination of ceramic diaphragm sensor and ceramic jacketed Pirani sensor ensures reliable readings in a wide range.*

# EXCLUSIVE DEAL - SAVE OVER £300